


Vascular Intervention // **Peripheral**
Drug-Coated Balloon Catheter/0.018"/OTW

Passeo[®]-18 Lux[™]



Clinically proven



Effective drug delivery



Prolonged drug presence



Passeo-18 Lux

Indicated to dilate de novo or restenotic lesions in the infrainguinal arteries.*

Technical Data	Drug-coated balloon
Catheter type	OTW
Recommended guide wire	0.018"
Tip	Short, tapered
Balloon markers	2 swaged markers (zero profile)
Shaft	3.8F, hydrophobic coated
Usable Length	90, 130 cm; 150 cm (only \varnothing 2.0 mm)
Introducer size	4F (\varnothing 2.0 - 4.0 mm); 5F (\varnothing 5.0 - 7.0 mm)
Nominal Pressure (NP)	6 atm
Rated Burst Pressure (RBP)	15 atm (\varnothing 2.0 - 5.0 mm); 12 atm (\varnothing 6.0 - 7.0 mm)

Coating

Drug	Paclitaxel
Drug concentration	3.0 $\mu\text{g}/\text{mm}^2$
Coating matrix	Paclitaxel and Butyryl-tri-hexyl citrate (BTHC)
Coated area	Cylindrical section of the balloon, exceeding the proximal and distal markers

Compliance Chart

Balloon diameter x length (mm)

		\varnothing 2.0 x 40-120	\varnothing 2.5 x 40-120	\varnothing 3.0 x 40-120	\varnothing 4.0 x 40-120	\varnothing 5.0 x 40-120	\varnothing 6.0 x 40-120	\varnothing 7.0 x 40-120
Nominal Pressure (NP)	atm**	6	6	6	6	6	6	6
	\varnothing (mm)	2.0	2.5	3.0	4.0	5.0	6.0	7.0
Rated Burst Pressure (RBP)	atm**	15	15	15	15	15	12	12
	\varnothing (mm)	2.1	2.6	3.2	4.3	5.3	6.2	7.3

**1 atm = 1.013 bar

Ordering Information

Catheter Length (cm) Balloon \varnothing (mm) Balloon Length (mm)

		Balloon Length (mm)			
		40	80	120	
4F	90	2.0	379860	379861	379862
	90	2.5	379866	379867	379868
	90	3.0	370843	370848	370853
	90	4.0	370844	370849	370854
	90	5.0	370845	370850	370855
5F	90	6.0	370846	370851	370856
	90	7.0	370847	370852	370857
	150	2.0	379863	379864	379865
4F	130	2.5	379869	379870	379871
	130	3.0	370858	370863	370868
	130	4.0	370859	370864	370869
	130	5.0	370860	370865	370870
5F	130	6.0	370861	370866	370871
	130	7.0	370862	370867	370872

*Indication as per IFU.

Passeo and Lux are trademarks or registered trademarks of the BIOTRONIK Group of Companies.